



NexxZr® Shine

Effective polishing.  
Natural gloss.

Effective polishing.  
Natural gloss.

# NexxZr® Shine



**NexxZr® Shine is a two-step diamond polishing system for the effective polishing of dental ceramics, e.g. zirconium oxide, aluminium oxide, lithium disilicate or glass ceramic.**

The dimensionally stable, diamond-infused rubber polishers convince users with their easy and efficient application. Additionally, they have a long service life and provide a true-to-nature lustre. NexxZr Shine creates a homogeneous

surface with both monolithic and partially veneered restorations. There are three coordinated shapes for the optimum polishing of the different structures: Turn 18, Twist 17 and Tip 5. The shapes are available in the grits smooth (Step 1, prep polishing) and gloss (Step 2, high-gloss polishing). The universally usable diamond polishing paste NexxZr Shine Paste completes the assortment. It removes micro-scratches and leaves a smooth-as-glass surface with a high gloss.

## Technical data



**Type**  
Two-step diamond  
polishing system



**Material**  
Diamond  
infused rubber  
polishers



**Grits**  
smooth, gloss



**Suitable for**  
all ceramics



**3 molds, 1 paste**  
Turn 18, Twist 17,  
Tip 5, Paste

### 1 Two-step diamond polishing system

Easy and efficient application on all ceramic materials




### 2 Dimensionally stable, diamond-infused rubber polishers

True-to-nature lustre with a long service life



### 3 Three coordinated shapes and a paste

Optimum polishing results, irrespective of structure and surface



# Delivery forms



Kit

**NexxZr Shine Starter Kit**  
Turn 18 smooth / 1pcs.  
Turn 18 gloss / 1pcs.  
Twist 17 smooth / 2pcs.  
Twist 17 gloss / 2pcs.  
Tip 5 smooth / 2pcs.  
Tip 5 gloss / 2pcs.  
Paste 20 g / 1pcs.

**NexxZr Shine Polishing Kit**  
Turn 18 smooth / 1pcs.  
Turn 18 gloss / 1pcs.  
Twist 17 smooth / 2pcs.  
Twist 17 gloss / 2pcs.  
Tip 5 smooth / 2pcs.  
Tip 5 gloss / 2pcs.



Refill

NexxZr Shine Turn 18 smooth / 5pcs.  
NexxZr Shine Turn 18 gloss / 5pcs.  
NexxZr Shine Twist 17 smooth / 5pcs.  
NexxZr Shine Twist 17 gloss / 5pcs.  
NexxZr Shine Tip 5 smooth / 5pcs.  
NexxZr Shine Tip 5 gloss / 5pcs.  
NexxZr Shine Paste 20 g / 1pcs.

## Manufacturer

**Sagemax Bioceramics, Inc.**  
34210 9th Ave. South, Suite 118  
Federal Way, WA 98003, USA  
[sagemax.com](http://sagemax.com)