



NexxZr.T Multi & NexxZr.+ Multi 2.0

Explore multi zirconia.  
Experience esthetics.

Explore multi zirconia.  
Experience esthetics.

# NexxZr® T Multi



Esthetics



Range of indications



1

## Smooth color and translucency gradient

Multifunctional esthetics

2


## High strength

Wide variety of indications

3

## High-quality production process

Tension-free sintering without distortion  
and high accuracy of fit



**Type / Class:**

Type II / Class 5

**Indication:**

- › Single-unit restorations
- › 3-unit bridges
- › Multi-unit bridges with max. 2 connected pontics

**Processing techniques:**

- › Infiltration technique
- › Staining technique
- › Cut-back technique

**Geometries:**

W-98, Z-95

**Thicknesses:**

16, 20, 25 mm

**Colors:**

BL1, BL2, BL3, BL4, A1, A2, A3, A3.5, A4,  
B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4



# NexxZr<sup>®</sup>+ Multi 2.0

Esthetics

\*\*\*

Range of indications

\*\*

1

## Natural color gradient

Efficient fabrication of natural looking restorations

2

## Specially adjusted translucency gradient

Good coverage in the cervical area and increased translucency in the incisal area

3


## High flexural strength

Single-unit restorations up to 3-unit bridges

5Y-TZP  
650 MPa  
49% Translucency

Average values

4Y-TZP  
850 MPa  
45% Translucency



Incisal  
3 mm

Mixture  
4 mm

Dentin  
7, 9 or 13 mm  
depending on  
disc thickness

**Type / Class:** Typ II / Class 4

**Indication:**

- › Single-unit restorations
- › 3-unit bridges

**Processing techniques:**

- › Infiltration technique
- › Staining technique
- › Cut-back technique

**Geometries:** W-98


**Thicknesses:** 14, 16, 20 mm

**Colors:** BL1, BL2, BL3, BL4, A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4

## NexxZr T Multi



Single-unit restoration up to multi-unit bridges - all with one material



### NexxZr T Multi

A multilayer zirconium oxide with multifunctional esthetics and a wide variety of indications



## NexxZr+ Multi 2.0



Natural color and translucency gradient



### NexxZr+ Multi 2.0

Multilayer zirconium oxide with natural esthetics



### Manufacturer

**Sagemax Bioceramics, Inc.**  
34210 9th Ave. South, Suite 118  
Federal Way, WA 98003, USA  
[sagemax.com](http://sagemax.com)