



Paint 3D

Universal coloring system.
3D effect pastes.

Universal coloring system.
3D effect pastes.

Paint 3D



Paint 3D is a universal coloring system for the individual staining and glazing of dental ceramics in the CTE range from 9.4 to $17.5 \times 10^{-6}/K$.

Ideally suited for milling, pressing and veneering ceramics – including zirconium oxide, lithium disilicate and other glass ceramics – Paint 3D offers a wide range of applications and maximum flexibility in everyday laboratory work.

The innovative 3D effect pastes open up new creative possibilities and give surfaces an impressive three-dimensional depth and texture. Whether for subtle accents or striking effects: These pastes make it possible to transform even the smallest details into realistic dental works of art.

Thanks to their special consistency, the pastes can be applied and modeled effortlessly, making them ideal for demanding projects where precision and efficiency are paramount. Paint 3D supports time-saving workflows without compromising on esthetics or attention to detail – for professional results at the highest level.

Technical data



Typ
Low-fusing glass
ceramic



Pastes
Self-glazing
3D effect pastes



CTE range
 $9.4 - 17.5 \times 10^{-6}/K$



Applications
Zirconium oxide,
Lithium-Disilicate,
Veneering ceramics



Application techniques
Wet-in-Wet,
Classic


Highlights

1 Universal coloring system

Versatile application on dental ceramics in the CTE range of 9.4 to $17.5 \times 10^{-6}/K$



Applicable to zirconia, milled ceramics, pressed ceramics, and layered ceramics (full-ceramic & metal-ceramic)



Firing-stable colors – what you see is what you get!

2 3D effect pastes

Creates three-dimensional effects and textures for a lively, natural esthetics



Unique optical depth, even at minimal layer thickness



Special texture paste for voluminous, three-dimensional surface textures

3 Easy handling

Intuitive and effortless application for time-saving and consistent results



Ready-to-use consistency



Self-glazing properties for easy application

Delivery forms

Colors
White Esthetics - pastes with fluorescence

Dentin

767480	767481	767482	767483	767484

Paste per jar
4g

Enamel

767485	767496	767497	767498	767499	767500

Effect

767501	767502	767503

Intense Effect

767504	767505	767506	767507	767508

Glaze & Texture

767510	767509

Pink Esthetics - pastes without fluorescence

Base

767513	767514	767515

Effect

767516	767517	767518	767519	767520	767521

Glaze & Texture

767523	767522

Manufacturer

DSSM AG
Im alten Riet 9
9494 Schaan
Liechtenstein
sagemax.com

Universal Liquid

Starter Kit

Gingiva Kit

767512	767479	767512