

sagemax.com

sagemax®



Sagemax experience

---

NexxZr+ Multi 2.0  
micro cut-back  
by Mikhail Leushin

NexxZr+ Multi 2.0 micro cut-back  
by Mikhail Leushin

---

# Sagemax experience

## Initial situation

---







The initial situation was highly unsatisfactory for the patient.



Her greatest wish was to finally have a radiant smile again.



## Design process

---



After an in-depth consultation with the patient, a comprehensive reconstruction of the upper and lower jaw was planned. Four implants were placed in each arch. The highly translucent zirconia NexxZr+ Multi 2.0 was used as the framework material. The images show the CAD design of the framework.

## Framework finish



The zirconia restoration milled from NexxZr+ Multi 2.0 was minimally adjusted before the sintering process. The labial surfaces of the upper anterior teeth were slightly reduced and subsequently rebuilt after sintering using a micro-layering technique with IPS e.max Ceram materials (Co. Ivoclar). The gingiva-colored areas were also reconstructed with veneering ceramics.



The clinical final result demonstrates an outstandingly natural and highly esthetic outcome.



## Final result

---



Before.



After.



The patient was overjoyed with her new smile. She gained not only a new smile, but also a completely new sense of life. It's exactly these moments that make me love my job.



### Dental designer

**Mikhail Leushin**  
Dental Technician, Dubai

📷 skillslabdubai  
📱 Skills Tech Dental Lab  
[www.skillsdl.com](http://www.skillsdl.com)

---

**Contact us**  
[info@sagemax.com](mailto:info@sagemax.com)

[sagemax.com](http://sagemax.com)

**Follow us**  
[@sagemax.official](https://www.instagram.com/sagemax.official)  
📺 📺 📺 📺