

sagemax.com

sagemax®



Sagemax experience

---

NexxZr<sup>+</sup> Multi  
stained and micro-layered  
by Dominik Monreal

NexxZr<sup>+</sup> Multi stained and micro-layered  
by Dominik Monreal


# Sagemax experience

Initial situation



Heavy abrasion and deformation of the teeth. Patient's request: Complete restoration of the upper and lower jaw.


Framework fabrication, staining and micro-layering



Zirconium oxid: NexxZr<sup>+</sup> Multi.



Separated restorations.



Restorations before sintering.



Partially reduced anterior crowns and monolithic posterior bridges after sintering.



Monolithic restorations of the lower jaw after sintering.



Partially reduced anterior crowns.



Wash firing with IPS e.max Ceram, Co. Ivoclar.




Minimal dentine and incisal build-up with IPS e.max Ceram, Co. Ivoclar.



Elaboration of the surface structure.



Staining and glazing with IPS Ivocolor, Co. Ivoclar.



Stain and glaze firing.



Finished upper jaw restorations.



Finished lower jaw restorations.



Lifelike looking lower jaw restorations, ready for cementation.

## Final result



Before.



After.



Excellent final esthetic result and an extremely happy patient.

### Dental designers



**Dominik Monreal**  
Dental Technician, Germany

📷 dominikmonreal8



**Kristina Günter**  
Dentist, Germany

www.zahnarzt-kleve.de  
📘 Praxis für Zahnheilkunde Kristina Günter  
📷 zahnarzt\_kleve

Contact us  
info@sagemax.com

sagemax.com

Follow us  
@sagemax.official

