

sagemax.com

sagemax®



Sagemax experience

---

NexxZr Shine  
polishing  
by Ivan Unfeter

NexxZr Shine polishing  
by Ivan Unfeter

# Sagemax experience

Initial situation



Heavily discolored and deformed teeth. Patient's request: A completely new, beautiful smile.

Milling, texturing & sintering



Manual elaboration of the macro- and microtexture in the raw zirconium oxide restorations.



## Polishing, staining & glazing

---



Crowns with preserved texture and surface after sintering.



Polishing of all concave and convex surfaces with NexxZr Shine polishers. The texture remains present and visible.



Characterization and glazing with IPS Ivocolor (Co. Ivoclar).



The same process is carried out for the lower jaw restorations.

## Final result

---



Before.



After.





The final result fits perfectly into the surrounding oral cavity and gives the patient a new smile.



### Dental designer

**Ivan Unfeter**  
Dental Technician, Serbia

 dent.artist.unfeter  
 dent.artist.unfeter

---

**Contact us**  
info@sagemax.com

sagemax.com

**Follow us**  
@sagemax.official

