

sagemax.com

sagemax®



Sagemax experience

NexxZr Shine
Set za poliranje
by Ivan Unfeter

NexxZr Shine Set za poliranje
by Ivan Unfeter

Sagemax experience

Prvobitna situacija



Diskoloracije i narušena anatomija i izgled zuba. Zahtev pacijenta: Novi, lepši i blistavi osmeh.

Rezanje, obrada i sinterovanje



Manuelna obrada i definicija makro i mikro teksture, na sirovim Zirconium-oxide krunama, odmah posle rezanja.

Poliranje, karakterizacija i glaziranje




Krone sa površinskom teksturom nakon sinterovanja.



Poliranje konveksnih i konkavnih površina sa NexxZr Shine instrumentima za poliranje. Tekstura ostaje prisutna i vidljiva.



Karakterizacija i glaziranje sa IPS Ivocolor (Co. Ivoclar).



Isti proces je primenjen i na kronicama donje vilice.

Finalni rezultat



Pre.



Posle.





Finalni rezultat - krunice se savršeno uklapaju sa mekim tkivima i daju pacijentu novi blistavi osmeh.



Dental designer

Ivan Unfeter
Dental Technician, Serbia

 dent.artist.unfeter
 dent.artist.unfeter

Contact us
info@sagemax.com

sagemax.com

Follow us
@sagemax.official

