

sagemax.com

sagemax®



Sagemax experience

NexxZr T Multi
micro-layered
by Bastian Wagner

NexxZr T Multi micro-layered
by Bastian Wagner

Sagemax experience

Initial situation



Heavily discolored and deformed anterior teeth 11 and 21.



Patient's request: Masking of anterior teeth 11 and 21 with crowns.

Color determination with eLAB (Co. Emulation S.Hein)




Reflective image.



Polarized image.

Micro-layering with IPS e.max Ceram and IPS Ivocolor (Co. Ivoclar)



Design of the micro-layering. NexxZr T Multi A2 framework layer thickness approx. 0.6 mm.

Cementation



Micro-layered NexxZr T Multi crown during cementation.

Excellent coverage of the dark tooth stump.

Final result



Before.



After.




The final result fits perfectly into the surrounding oral cavity.



Dental designer

Bastian Wagner
Master Dental Technician, Germany

www.wagner-dental-design.com

 [bastian.wagner.12](https://www.facebook.com/bastian.wagner.12)

 [mdt_bastian_](https://www.instagram.com/mdt_bastian_)

Contact us
info@sagemax.com

sagemax.com

Follow us
[@sagemax.official](https://www.instagram.com/sagemax.official)

