

sagemax.com

sagemax[®]






Sagemax experience

NexxZr T Multi
monolithic
by Mateja Stevanović

NexxZr T Multi monolithic
by Mateja Stevanović


Sagemax experience

Initial situation



Dissatisfied with the appearance of her teeth, the patient visited our practice with the specific desire to have the existing diastemas closed. At the same time, she emphasized the importance of preserving the natural shape and overall aesthetic of her smile.

Digital design & material selection



After completing the 3D smile design, the final design is ready for milling.




NexxZr T Multi A1 from the Sagemax portfolio was chosen as the ideal solution due to its high translucency and optimized strength – 1200 MPa in the cervical area and 650 MPa in the incisal area – which helps prevent wear on the antagonists. The crowns were milled in a full-anatomy design.

Fabrication of the restorations



In this phase, we manually refine the details, adding texture and morphology, after which the crowns are ready for the sintering process.



During the characterisation phase, we bring the monolithic crowns to life, adding depth and vitality. Here, we used pigments from the Ivocolor Essence palette, then finished with Ivocolor Fluo Glaze (Co. Ivoclar).



Final result



The combination of all materials led to a stunning result. Even with a difficult tooth position, the outcome was a naturally beautiful and harmonious smile.



Dental designer

Mateja Stevanović
Dental Technician, Serbia

📷 [tehnodent.vranje](#)

Contact us
info@sagemax.com

sagemax.com

Follow us
[@sagemax.official](#)

