

sagemax.com

sagemax®



Sagemax experience

NexxZr T Multi
monolithic
by Michał Nytko



NexxZr T Multi monolithic
by Michał Nytko

Sagemax experience

Initial situation




Heavily discolored and deformed anterior teeth. Patient's request: masking of the anterior teeth with crowns.




Mock-up to define the midline and length of the teeth.


Framework design




A pre-op scan was used to select the appropriate shape of the new teeth.



Framework design of the restorations.




Framework design of the restorations.



Space conditions.

Fabrication of the restorations



The restorations were fabricated fully monolithically from NexxZr T Multi A1. The restorations were then stained with Creation Magic Color (Co. Creation Willi Geller) and glazed with NexxZr Glaze Spray.

Finished NexxZr T Multi restorations



Smooth color and translucency gradient with NexxZr T Multi.



Smooth, homogeneous surfaces with NexxZr Glaze Spray. Fully monolithic, only stained and glazed.

Final result



Before.



After.





The restorations give the patient a new smile.



Dental designer

Michał Nytko
Dental Technician, Poland

 dentnytko
 michael_nytko

Contact us
info@sagemax.com

sagemax.com

Follow us
@sagemax.official

