



# sagemax®

## Zirconia Portfolio

# NexxZr®



Product	NexxZr S	NexxZr T	NexxZr T Multi	NexxZr+	NexxZr+ Multi
Strength	1370 MPa	1270 MPa	Incisal 630 MPa Dentin 1170 MPa	Preshaded 880 MPa White 1000 MPa	Incisal 630 MPa Dentin 880 MPa
Translucency	30%	42%	Incisal 49% Dentin 42%	46%	Incisal 49% Dentin 46%
Colors	White	White, OM1, OM2, OM3, A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4	BL1, BL2, BL3, BL4, A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4	White, A1, A2, A3, A3.5, B1, B2	A1, A2, A3, B1, B2, C2, D2
Geometries	W-98, Z-95	W-98, Z-95	W-98, Z-95	W-98, Z-95	W-98, Z-95
Thicknesses	10, 12, 14, 16, 18, 20, 22, 25 mm	10, 12, 14, 16, 18, 20, 22, 25, 30 mm	16, 20 mm	10, 12, 14, 16, 18, 20, 22, 25, 30 mm	16, 20 mm

sagemax.com